



- DPX-170S features include;
- Variable Speed Control offers maximum precision & accuracy
 - High Powered pneumatic motor
 - Fully Autoclavable
 - Low Maintenance



Pneumatic Saw Specification

Saw model	Order code	Speed	Oscillating arc	Operating pressure
DPX-170S	224774	0 - 17000 cpm	6°	6 - 7 bar

Wash Assist Hose

Model	Order code	Length	Probe
Air hose	13000	3m	MA7 Schrader
Air hose	9500	3m	1/4 NTP Connection

Other probe connectors are available on request



A range of TPLO saw blade sizes are available for purchase separately.

* Drive adaptor Slocum 16330 and 4mm hex key 670423 sold separately